

QTL Analysis Pipeline - Use Cases Guide

Quick Start Examples

Basic eQTL Analysis

Use Case: Standard expression QTL analysis with VCF genotypes

bash

- *# 1. Edit main config with your file paths*
- `nano config/config.yaml`
-
- *# 2. Update these paths in config.yaml:*
- `input_files:`
- `genotypes: "your_data/genotypes.vcf.gz"`
- `covariates: "your_data/covariates.txt"`
- `annotations: "your_data/annotations.bed"`
- `expression: "your_data/expression.txt"`
-
- *# 3. Run eQTL analysis*

```
python run_QTLPipeline.py --config config/config.yaml --analysis-types eqtl
```

Expected Output: `results/qlt_results/eqtl_significant.txt` with significant eQTL associations

Use Case 1: GATK4.5 Joint Calling VCF

Scenario

You have genotype data from GATK4.5 joint calling in VCF.GZ format with sample names like `sample_001_L001`.

Configuration

bash

- *# Use the GATK-optimized configuration*
- `cp config/gatk_joint_calling.yaml config/my_gatk_analysis.yaml`
-
- *# Edit with your paths*

`nano config/my_gatk_analysis.yaml`

Key Settings for GATK:

yaml

- `genotype_processing:`
- `normalize_chromosomes: true`
- `chromosome_prefix: "chr"` *# Converts "1" → "chr1"*
- `handle_multiallelic: true` *# Splits multi-allelic sites*
- `remove_phasing: true` *# Removes phasing information*






`min_call_rate: 0.98` *# GATK has high call rates*

Run Command

bash

`python run_QTLPipeline.py --config config/my_gatk_analysis.yaml`

Expected Processing

-  Automatic format detection
 -  Chromosome naming normalization
 -  Multi-allelic site splitting
 -  Phasing removal
 -  Quality filtering
-

Use Case 2: PLINK BED Format

Scenario

You have genotype data in PLINK BED format (.bed, .bim, .fam files).

Configuration

bash

- `cp config/plink_format.yaml config/my_plink_analysis.yaml`

`nano config/my_plink_analysis.yaml`

Key Settings for PLINK:

yaml

- `input_files:`
- `genotypes: "your_data/plink_data.bed" # .bim and .fam must exist`
- `genotype_processing:`
- `auto_detect_format: true # Automatically converts to VCF`




`chromosome_prefix: "none" # PLINK uses "1" not "chr1"`


Run Command

bash

`python run_QTLPipeline.py --config config/my_plink_analysis.yaml`

Expected Processing

-  Automatic BED → VCF conversion
-  Sample information from .fam file
-  Variant information from .bim file

-  Standardized output format
-

Use Case 3: Multi-omics Analysis

Scenario

You want to run eQTL, pQTL, and sQTL analyses together to compare molecular QTLs.

Configuration

bash

- `cp config/multi_omics.yaml config/my_multiomics.yaml`

`nano config/my_multiomics.yaml`

Key Settings for Multi-omics:

yaml

- `input_files:`
- `expression: "data/expression.txt"`
- `protein: "data/protein.txt"`
- `splicing: "data/splicing.txt"`
-
- `analysis:`
- `qtl_types: "all"` *# Runs all three analyses*
-
- `plotting:`
- `plot_types:`
- - `"manhattan"`
 - `"qq"`
-
- `"summary"` *# Cross-analysis comparison plot*

Run Command

```
bash
```

```
python run_QTLPipeline.py --config config/my_multiomics.yaml
```

Expected Outputs

- `eql_significant.txt` - Expression QTLs
 - `pqtl_significant.txt` - Protein QTLs
 - `sqtl_significant.txt` - Splicing QTLs
 - `analysis_summary.png` - Comparison of all three
-

Use Case 4: GWAS + QTL Integration

Scenario

You want to run both QTL analysis and GWAS on the same genotype data.

Configuration

```
bash
```

- `cp config/multi_omics.yaml config/my_gwas_integration.yaml`

```
nano config/my_gwas_integration.yaml
```

Key Settings for GWAS:

```
yaml
```

- `analysis:`
- `qtl_types: "eqtl"` *# Run eQTL analysis*
- `run_gwas: true` *# Enable GWAS*
- `gwas_phenotype: "data/gwas_phenotypes.txt"`
-
- `gwas:`
- `method: "linear"` *# For continuous traits*

- `# method: "logistic"` `# For case-control traits`

```
covariates: true # Adjust for covariates
```

GWAS Phenotype File Format:

text

- | sample_id | height | weight | disease_status |
|-----------|--------|--------|----------------|
| sample1 | 175.2 | 68.5 | 0 |
| sample2 | 168.7 | 72.1 | 1 |
-
-

| | | | |
|---------|-------|------|---|
| sample3 | 182.3 | 75.8 | 0 |
|---------|-------|------|---|

Run Command

bash

```
python run_QTLPipeline.py --config config/my_gwas_integration.yaml
```

Expected Outputs

- QTL results in `qtl_results/`
 - GWAS results in `gwas_results/`
 - Combined plots showing both analyses
-

Use Case 5: Large-Scale Study

Scenario

You have thousands of samples and need to optimize for computational efficiency.

Configuration

bash

- `cp config/large_scale_eqtl.yaml config/my_large_study.yaml`

`nano config/my_large_study.yaml`

Optimization Settings:

yaml

- `qtl:`
- `permutations: 1000` *# Reduced for efficiency (default: 1000)*
- `# permutations: 100` *# For very large studies*
-
- `output:`
- `remove_intermediate: true` *# Save disk space*
- `compression: true` *# Compress large files*
-
- `genotype_processing:`
- `filter_variants: true`
- `min_maf: 0.01` *# Filter rare variants to reduce load*

`min_call_rate: 0.95`

Run Command

bash

- *# Test with subset first*
- `python run_QTLPipeline.py --config config/my_large_study.yaml`
`--analysis-types eqtl`
-
- *# Full run*

`python run_QTLPipeline.py --config config/my_large_study.yaml`

Performance Tips

- Start with fewer permutations for testing
- Use `--validate-only` first to check inputs

- Monitor memory usage in large runs
 - Consider running on HPC cluster for very large datasets
-

Use Case 6: Fine-Mapping Analysis

Scenario

You want high-resolution QTL mapping with strict quality controls.

Configuration

bash

- `cp config/fine_mapping.yaml config/my_finemapping.yaml`

`nano config/my_finemapping.yaml`

Fine-Mapping Settings:

yaml

- `qtl:`
- `cis_window: 500000` *# Smaller window for fine-mapping*
- `permutations: 10000` *# More permutations for accuracy*
- `fdr_threshold: 0.01` *# Stricter significance threshold*
- `genotype_processing:`
- `min_maf: 0.01` *# Include rare variants*
- `min_call_rate: 0.99` *# High quality requirement*

`quality_threshold: 50` *# Strict quality filter*

Run Command

bash


```
python run_QTLPipeline.py --config config/my_finemapping.yaml
```

Expected Results

- Higher precision QTL mapping
 - More accurate FDR estimates
 - Better variant-gene pairing
-

Use Case 7: Tissue-Specific Analysis

Scenario

You have data from specific tissues and want tissue-specific QTLs.

Configuration

yaml

- *# Start with main config and modify:*
- **qtl:**
- **cis_window:** 500000 *# Smaller window for tissue-specific*
-
- **genotype_processing:**
- **filter_variants:** true
- **min_maf:** 0.05 *# Higher MAF for power in smaller samples*
-
- **analysis:**

```
qtl_types: "eqtl" # Focus on expression
```

Run Command

bash

```
python run_QTLPipeline.py --config config/config.yaml --analysis-types eqtl
```

Troubleshooting Common Scenarios

Sample Name Mismatch

Problem: Samples in genotype file don't match phenotype file

Solution:

yaml

- `genotype_processing:`
- `match_samples: true`

```
sample_matching_method: "prefix" # or "pattern"
```

Memory Issues

Problem: Pipeline runs out of memory

Solution:

yaml

- `qtl:`
- `permutations: 100` *# Reduce permutations*
- `output:`

```
remove_intermediate: true # Clean up files
```

Long Runtime

Problem: Analysis taking too long

Solution:

yaml

- `qtl:`
- `permutations: 1000` *# Reduce from default*
- `cis_window: 500000` *# Smaller window*
- `genotype_processing:`

`min_maf: 0.05` *# Filter more variants*

Quick Reference Commands

| Use Case | Command |
|-----------------|---|
| Basic eQTL | <code>python run_QTLPipeline.py --config config.yaml --analysis-types eqtl</code> |
| GATK VCF | <code>python run_QTLPipeline.py --config gatk_joint_calling.yaml</code> |
| PLINK BED | <code>python run_QTLPipeline.py --config plink_format.yaml</code> |
| Multi-omic s | <code>python run_QTLPipeline.py --config multi_omics.yaml</code> |

| | |
|--------------|--|
| With GWAS | <code>python run_QTLPipeline.py --config config.yaml --run-gwas</code> |
|--------------|--|

| | |
|------------------|---|
| Validate Only | <code>python run_QTLPipeline.py --config config.yaml --validate-only</code> |
|------------------|---|

| | |
|----------------|---|
| Large Study | <code>python run_QTLPipeline.py --config large_scale_eqtl.yaml</code> |
|----------------|---|

-